

HANSATON wave

**Help your clients
embrace life's beauty**

Your clients can enjoy the benefits
of connectivity to their
mobile phones along with
the sounds they love.



Meeting the needs of your clients

Boosts conversations and suppresses noise

HearIntelligence™ technology supports your clients' hearing by processing incoming sounds and reducing them to essential components that are either boosted or suppressed, providing an optimized hearing experience.


Design meets performance

HANSATON combines high performance with stylish design in each of the HANSATON wave styles. You can select from a RIC design that is very high-fashion or BTEs that are known for their reliability. Stylish designs can make your clients feel comfortable trying hearing aids.

Sound Optimization

With the FocussedFit workflow you can get your clients on their way to hearing better quickly. The AcclimatizationManager starts them at a level that is comfortable and provides a way to reach an optimal setting over time as they become accustomed to amplification. The AutomaticMic is an easy solution, which automatically changes directionality based on the noise level, so your clients can enjoy their hearing aids automatically.

Encourage your clients to hear the important moments in their lives.



Wireless connectivity features

Hands-free phone calls and streaming with all mobile phones or other Bluetooth® devices.*

Easy connectivity options to share with your clients:




Connections to all mobile phones

Turn their hearing aids into discrete wireless headphones. Enable your clients to answer a phone call or listen to music streaming directly to their hearing aids and to connect two phones or tablets at once.*

Personalized hearing experience

Your clients can make their own personalization choices in the moment with our smartphone App, which has a nice balance of quick meaningful adjustments and more advanced fine-tuning options.



TV Connector

Clients can watch their favorite movies, TV programs or catch the news. The TV Connector sends the audio from their TV or other audio devices directly to their hearing aids.

*Smartphones and traditional mobile phones with compatible Bluetooth® profiles.
The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.

Beautifully designed to meet the diverse needs of your clients

HANSATON wave offers 1 RIC and 3 BTE hearing aid solutions. This enables you to fit a range of hearing losses – from mild to profound – to meet the needs of more clients.

Based on proven technology, HANSATON wave covers a wide range of hearing losses and has appealing designs to meet your client's needs.



Water and dust resistant





Equipped with a special plasma coating, these hearing aids are elegant in appearance and are dust and water resistant according to the stringent IP68 classification.

HANSATON wave R312

RIC hearing aid

The HANSATON wave R312 hearing system is for mild to profound hearing losses.




	S Receiver	M Receiver	P Receiver	UP Receiver
Output / Gain (2cc)	111 / 47	114 / 51	122 / 59	132 / 71
Fitting ranges				

HANSATON wave 312 M

BTE hearing aid

The HANSATON wave 312 M hearing system is for mild to severe hearing losses.



	M Receiver
Output / Gain (2cc) Earhook	128 / 63
Output / Gain (2cc) Slim tube	122 / 58
Fitting range	

HANSATON wave 13 SP BTE hearing aid

The HANSATON wave 13 SP hearing system is for mild to severe hearing losses.



Battery Size 13
Rocker
Integrated Telecoil
Pre-sets in scout

SP Receiver	
Output / Gain (2cc) Filtered Earhook	133 / 75
Output / Gain (2cc) Unfiltered Earhook	139 / 81
Fitting range	

HANSATON wave 675 UP BTE hearing aid

The HANSATON wave 675 UP hearing system is for moderate to profound hearing losses.



Battery Size 675
Rocker
Integrated Telecoil
Pre-sets in scout

UP Receiver	
Output / Gain (2cc) Filtered Earhook	134 / 78
Output / Gain (2cc) Unfiltered Earhook	141 / 84
Fitting range	

HEAR Intelligence™

HearIntelligence™ technology was developed by HANSATON to support natural hearing by enhancing natural functions and processes in human hearing.

Focusing on Speech

- ✓ AutomaticMic
- ✓ ConversationLift
- ✓ SoundRestore frequency compression

Comfort in Noise

- ✓ NoiseReduction
- ✓ Sound Impulse Manager
- ✓ FeedbackManager

Personalized Adjustments

- ✓ 12 Channels
- ✓ FocussedFit
- ✓ Intelligent Acclimatization
- ✓ DataLogging
- ✓ Up to 4 manual programs

Wireless Connectivity

- ✓ Made for all phones
- ✓ Handsfree Bluetooth® phone calls
- ✓ Bluetooth media streaming
 - SpeechStream
 - MusicStream
- ✓ Pair up to 8 Bluetooth devices
- ✓ Connect up to 2 Bluetooth devices at once
- ✓ Optional Accessories:
 - TV Connector
 - PartnerMic
 - Remote Control
- ✓ stream remote App compatible
- ✓ BiLink
- ✓ Binaural connection loss indicator





HANSATON is the exclusive brand for you and your clients since 1957.

We focus on personal relationships and your success.

Our German heritage, combined with the innovative power of Sonova, enables us to create beautifully designed products, based on cutting-edge, proven technology.



We drive the transition to a resource-efficient and low-carbon future of our operations throughout the lifecycle of our products and services.

sonova.com/sustainability