

Easy Line ChargerGo ITE

Technical data

Dimensions (L x W x H)	99 x 57 x 54 mm (3.9 x 2.2 x 2.1") +/- 5%
Weight	84 g (3 oz) +/- 10%
Color	Pewter
Operating conditions	+5° to +40° Celsius (+41° to +104° Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during transportation and storage	0% to 93% (non-condensing)
Humidity during operation	0% to 90% (non-condensing)
Atmospheric pressure during use	500 hPa to 1060 hPa
Atmospheric pressure during storage and transportation	500 hPa to 1060 hPa
Included parts	Power supply USB-C cable
External interfaces	USB-C port
Battery	Rechargeable lithium-ion
Charging time	See user guide
Input	5 V DC +/-5%, current 500 mA
Rechargeable lithium-ion battery capacity	1.184 Wh

The battery cannot be overcharged. It may become warm during charging.

